



West Nile Virus (WNV) detection & quantitation by qRT-PCR (Eva green).

Created: Oct 18, 2016; Last modified: Aug 11, 2022; Version: 4.0

This protocol describes the real-time polymerase chain reaction with reverse transcription (qRT-PCR) approach developed for the detection and quantitation of viral titres in biological specimens using either a two-step (separate RT and PCR reactions) or one-step EVA-Green modality (but easily TaqMan adaptable). WNV is a group IV (+)ssRNA species of an unassigned order, Flaviviridae family, Flavivirus genus related to the dengue, yellow fever, Japanese encephalitis, and Zika viruses. WNV is a member of the Japanese encephalitis antigenic serocomplex of viruses, exhibits a similar virion morphology as DENV of 45-50 nm. Its genome is 11 to 12 kb in length encodes seven nonstructural proteins and three structural proteins. The RNA strand is held within a nucleocapsid formed from 12-kDa protein blocks; the capsid is contained within a host-derived membrane altered by two viral glycoproteins.

Oligonucleotide primers

Name	Sequence ¹	bp	%GC	Tm	Hair	HomD	HetD	Amplicon
WNV-F	5'-CAg-ACC-ACg-CTA-Cgg-Cg-3'	17	70.6	59	-3.1	-6.75	-13.12	101
WNV-R	5'-CTA-ggg-CCg-CgT-ggg-3'	15	80	59.4	-0.59	-10.36		

Note 1: Reverse oligonucleotide primer sequences given in this table are the reverse-complement of sequence present in alignments and as they should be ordered for synthesis; Hairpin (Hair), homodimer (homD) and heterodimer (HetD) ΔQ 's given in kcal/mol.

Oligonucleotide primer map in artificial gene² (underlined primer sites)

AAGCTCCGCGTCTTTACCAAGGAAATAACATCACTGTAAGTGCCTATGCAAACGGCGACCATGCCGTACAGTTAAGGACGCCA
AATTCATTGTGGGGCCAATGTCTTCAGCCTGGACACCTTTCGACAACAAAATTTGTGGTGTACAAAAGGTGACGTCTATAATATGGA
CTACCCGCCCTTTGGCGCAGGAAGACCAGGACAATTTGG | CAGACCACGCTACGGCGTGCTACTCTGCGGAGAGTGCAGTCTGCG
ATAGTGCCCCAGGAGGACTGGGTAAACAAAGGCCAAACCAACGC CCCACGCGGCCCTAG | AAGGACTAGAGGTTAGAGGAGACCAT
GTCCATGTCAACCCACGGTCACTCATCTGGCATGTCTTGGGTCTATCCTGTCTGTCTCTACAGCATCATTCAGGCACAGAACGCC
| AAGGACTAGAGGTTAGAGGAGACCCCCCAACACAAAACAGCATATTGACGCTGGGAAAGACCAGAGATCCTGCTGTCTCTGC
AACATCAATCCAGGCACAGAGCGCC | AAATACACATACCAAAACAAAGTGGTGAAGGTTCTCAGACCAGCTGAAGGAGGAAAAC
AGTCATGGACATCATCTCAAGCAAGACCAGAGAGGGAGTGGA | ATGTATGTGAGTGCTGATGCCACGAAATGGTCACCAGGAGA
TAATTCGGCAAAGTTAAGAGATTCACACAGGCATTATATGATGGCTTGTGAGATGAGAA GTTAAAAATGTTGCGTTGTTGATGC |
ATGTAACACCTCTACAATGATGCCGACAAGATTGTATTCAAAGTCAATAATCAGGTGGTCTCTTTGAAGCCTGAGATTATCGTG
GATCAA TATGAGTACAAGTACCCTGC | GCATAGCTGTAGGAAGGACTCTCAA GGCATTTTTCTCTGCTGCTTACTGACCCCTTT
AGGGAACGAAGCCCTGGGGATACTGTCTTGAAGTGGATGCTT | GCAACGCGGATTGAGTTCTCTTCACATAATCGCCCC
GAGCTCGCTTATCGTTAAGCAGCTCTGCGCTACTATGGGTCCCGTGTAGAGGC | TGATGATGCCGTCGTGTGCTACAA CAGTAA
CTATGCGGCTCAAGGTTTAGTAGCTAGCATTAAGAAGCTTAAAGGCAGTCTGTATTATCAAGATAATGTGTTTCATGTCTGAGGCA
AAATGTTGGACTGAGACTGACCTTACTAAA GGACCTCACGAATTTTGCTACA

Note 2: Primer sites for different pathogen sequences shown in red, TaqMan probe sites shown in blue,





WNV specific primer binding sites shown in grey highlight. WNV artificial nucleocapside 3' UTR gene sequence derived from Bethesda reference strain partial sequence as available from GenBank accession number AY646354.

Two-step (separate RT and PCR) components and conditions

Two-step (separate RT and PCR) uses M-MULV RT Enzyme (Vivantis Cat. No.: ME2305) and iTaq qPCR Master Mix (Biorad Cat. No.: 1725150).

First strand synthesis (RT)

	cf	1x
dH ₂ O	---	2.5 µL
10 µM Forward oligo	1.125 µM	2.3 µL
10 µM Reverse oligo	1.125 µM	2.3 µL
10 mM dNTPs 10 mM	250 µM	0.5 µL
RNA	---	10 µL
		vf: 17.6 µl

Run RT-1 program in Axygen TC-1

Total time: 6 min		
95 °C	2 min	1 cycle
4 °C	2 min	

	1x cf	
RT Buffer	1x cf	2 µL
RT Enzyme 250 IU/µl	5 IU/µL	0.4 µL
		vf: 20 µl

Run RT-2 program in Axygen TC-1

Total time: 1:12 hrs		
38 °C	60 min	1 cycle
95 °C	5 min	
4 °C	5 min	

Polymerase Chain Reaction (PCR)

	cf	1x
dH ₂ O	---	3.7 µL
2x iTaq SYBR Green	1 x	5 µL
10 µM Forward oligo	150 nM	0.15 µL
10 µM Reverse oligo	150 nM	0.15 µL
Template	10 ng	1 µL
		vf: 10 µl

Run generic program in Applied Biosystems 7500

Total time: 2:10 hrs		
95 °C	3 min	40 cycles
94 °C	15 sec	
60 °C ³	60 sec	
95 °C	15 sec	
60 °C	20 sec	
Ramp ³	60 min	
95 °C	15 sec	

Note 3: Data acquisition.



One-step qRT-PCR components and conditions

The One-Step qRT-PCR procedure uses SCRIPT One-Step RT-qPCR EvaGreenMaster (Jena Bioscience Cat No. PCR-514S).

qRT-PCR



		1x
dH ₂ O	---	1.2 μL
10 μM Forward oligo	400 nM	0.4 μL
10 μM Reverse oligo	200 nM	0.2 μL
ROX	0.5 x	0.2 μL
SCRIPT 2x	1 x	5 μL
Template	< 1 μg	3 μL
		vf: 10 μl

Run Green Master program in Applied 7500

50 °C	10 min	40 cycles
95 °C	5 min	
94 °C	15 sec	
60 °C ³	60 sec	
95 °C	15 sec	
60 °C	20 sec	
Ramp ³	60 min	
95 °C	15 sec	

Note 3: Data acquisition.

Performance summary



Standard curve: $m = -3.41$, $Y = 36.708$, $R^2 = 0.995$

Target amplicon T_m: $87.16^\circ\text{C} \pm 0.17 \text{ SD}$ (86.99 to 87.32°C)

Limit of detection (LODet): $7.51 \times 10^{-1} \text{ cp}/\mu\text{L}$

Limit of discrimination (LODis): $7.51 \times 10^2 \text{ cp}/\mu\text{L}$





Titration curve preparation for quantitative analysis of viral titres

Using Artificial Gene (AG) sequence stock having all viral references (CHIKV, WNV, DENV1-3, DENV2-4, ZIKV, SNV, RABV, TCRV, MERS and SARS) prepare the following log titres (stock currently at 7.51×10^7 cp/ μ L, working dilution of AG prepared at 7.51×10^6 cp/ μ L).

		AG stock	dH ₂ O vol	Ct
Working	7.5×10^7 cp/ μ L	10 μ L	90 μ L	9.20 (± 0.07)
1 st	7.5×10^6 cp/ μ L	10 μ L	90 μ L	12.78 (± 0.02)
2 nd	7.5×10^5 cp/ μ L	10 μ L	90 μ L	16.37 (± 0.29)
3 rd	7.5×10^4 cp/ μ L	10 μ L	90 μ L	20.25 (± 0.20)
4 th	7.5×10^3 cp/ μ L	10 μ L	90 μ L	22.85 (± 0.63)
5 th	7.5×10^2 cp/ μ L	10 μ L	90 μ L	27.15 (± 0.30)
6 th	7.5×10^1 cp/ μ L	10 μ L	90 μ L	29.32 (± 0.24)

Add 90 μ L to each of the 6 PCR 0.2 mL tubes. Take 10 μ L of initial working stock (at 7.51×10^7 cp/ μ L) and dispense into 1st PCR tube, wash tip 30 times, cap, vortex for 10 seconds and spin down for 10 seconds. Retrieve 10 μ L from volumetric centre of PCR tube and dispense into 2nd PCR tube repeating the same procedure for further dilutions.

References

1. Naze, F., Le Roux, K., Schuffenecker, I., Zeller, H., Staikowsky, F., & Grivard, P. et al. (2009). Simultaneous detection and quantitation of Chikungunya, Dengue and West Nile viruses by multiplex RT-PCR assays and Dengue virus typing using High Resolution Melting. *Journal Of Virological Methods*, 162(1-2), 1-7.

Revision history

- 1.0 Original document.
- 2.0 Optimized protocol.
- 3.0 Changes to document format only.
- 4.0 Changes to reagent volumes for optimized yield.





Standard Operating Procedures (SOPs)
Laboratorio de Genómica Viral y Humana
Facultad de Medicina UASLP



Distributed through a Creative Commons Attribution (BY) license granting the licensee the right to copy, distribute, display and make derivative works based on this document, including commercial use, as long as they credit the author as "*Laboratorio de Genómica Viral y Humana, Facultad de Medicina UASLP*".