



Human immunodeficiency virus (HIV) proviral DNA detection through end-point PCR.

Created: Aug 28, 2008; Last modified: Mar 23, 2021, Version: 4.0

This protocol describes the optimized end-point PCR component and conditions for the detection of integrated proviral sequences of the human immunodeficiency virus type 1. This method relies on a nested PCR approach in which the product of an initial PCR is used as a template for a second PCR using different (nested) oligonucleotides. This approach is more sensitive than single-pass PCRs and as such requires greater care and discipline to avoid contamination throughout setup.

Oligonucleotides

| Name | PCR | Sequence* | bp | %GC | Tm | Position † | Amplicon | Ref. |
|--------|-----|----------------------------------|----|-----|----|------------|----------|------|
| HIV-FO | 1 | 5'-TAC-Agg-AgC-AgA-TgA-TAC-Ag-3' | 20 | 45 | 50 | 141-161 | 294 | 1 |
| HIV-RO | | 5'-CCT-ggC-TTT-AAT-TTT-ACT-gg-3' | 20 | 40 | 48 | 418-438 | | |
| HIV-FI | 2 | 5'-ggA-AAC-CAA-AAA-TgA-TAg-gg-3' | 20 | 40 | 48 | 221-241 | 130 | |
| HIV-RI | | 5'-ATT-ATg-TTg-ACA-ggT-gTA-gg-3' | 20 | 40 | 48 | 331-351 | | |

* Reverse oligonucleotide primer sequences given in this table are the reverse-complement of the sequence present in the alignments and as they should be ordered for synthesis.

† Primer binding sites given on table are based on HXB2 reference sequence.

Components

| First PCR | | 1 rx (µl) |
|-------------------------|------------|-----------|
| dH ₂ O | Cf | 5.19 |
| 10x Buffer PCR | 1X | 1.25 |
| MgCl ₂ 50 mM | 3.0 mM | 0.75 |
| 4x dNTPs 10 mM | 200 µM | 0.25 |
| 2x Oligos 10 µM | 800 nM | 1.00 |
| Taq 5 UI/µL | 0.02 UI/µL | 0.063 |
| gDNA | - | 4.00 |
| | Vf | 12.5 |



Run HIV1

| Second (nested) PCR | | 1 rx (µl) |
|-----------------------------|------------|-----------|
| dH ₂ O | Cf | 8.44 |
| 10x Buffer | 1X | 1.25 |
| MgCl ₂ 50 mM | 2.0mM | 0.5 |
| 4x dNTPs 10 mM | 200 µM | 0.25 |
| 2x Oligos 10 µM | 800 nM | 1 |
| Taq 5 UI/µL | 0.02 UI/µL | 0.063 |
| 1 st PCR product | - | 1.00 |
| | Vf | 12.5 |



Run HIV2





Conditions (Total estimated time: 3 hrs)

| Time: 1:32 hrs | | Denaturing | | Annealing | Extension | | |
|----------------|-------------|------------|--------|-----------|-----------|-------|-------|
| HIV1 | Temperature | 94° | 94° | 55° | 72° | 72° | 4° |
| | Time | 2 min | 15 seg | 15 seg | 15 seg | 2 min | 5 min |
| | x30 cycles | | | | | | |

| Time: 1:32 hrs | | Denaturing | | Annealing | Extension | | |
|----------------|-------------|------------|--------|-----------|-----------|-------|-------|
| HIV2 | Temperature | 94° | 94° | 50° | 72° | 72° | 4° |
| | Time | 2 min | 15 seg | 15 seg | 15 seg | 2 min | 5 min |
| | x30 cycles | | | | | | |

Notes

1. Clean workbench with 0.1% NaOCl 0.1% followed by 70% Ethanol before and after work.
2. Preparation of RT mastermix should only be performed in the RT-PCR room.
3. Preparation of PCR mastermix and addition of sample DNA should only be performed in the pre-PCR enclosure or area.
4. Addition of positive template DNA should be performed on instrument (post-PCR) area.
5. All mastermixes should be prepared on ice to prevent excess evaporation.
6. Vortex and spin all mastermixes before and after aliquoting to PCR tubes.
7. For the first PCR distribute 8.5 μ L of the mastermix to each tube and then add 4 μ L of genomic DNA (\pm 100 ng/ μ L).
8. For the second PCR distribute 11.5 μ L of mastermix to each tube and then add 1 μ L of undiluted 1st PCR product.

References

1. Albert, J. Simple, sensitive, and specific detection of Human Immunodeficiency Virus type 1 in clinical specimens by polymerase chain reaction with nested primers. *Journal of Clinical Microbiology* **28**, 1560 - 1564 (1990).





Revision history

- 1.0 Original document.
- 2.0 Optimized protocol.
- 3.0 Optimized protocol
- 4.0 Changes to document format only.

