



Biological safety cabinet (BSC) use and knowledge assessment form. Created: Mar 28, 2022; Last modified: Mar 28, 2022, Version: 1.0

Name:							
Current aff	iliation:						
Instructor's name:							
1st test:	(dd/mmm/yyyy)	Score:	Certified:	□ Yes □	ı No		
2nd test: _	(dd/mmm/yyyy)	Score:	Certified:	□ Yes □	No No		
3rd test:	(dd/mmm/yyyy)	Score:	Certified:	□ Yes □	No		
4th test:	(dd/mmm/yyyy)	Score:	Certified:	□ Yes □	No		
Has subjec	t received previous BSC	training?	1) □ No (0)	Date:(dd/n	nmm/yyyy)		
Correctly d	lescribes main features of	\Box Yes $(+1^{ea})$	□ No (-1 ^{ea})	Pts:			
Correctly describes class II type A1, A2, B1, B2?			\Box Yes $(+1^{ea})$	□ No (-1 ^{ea})	Pts:		
Wears appropriate BSL2 personal protective equipment 2?			□ Yes (+1)	□ No (-1)	Pts:		
Wears appropriate BSL3 personal protective equipment 2?			□ Yes (+1)	□ No (-1)	Pts:		
Wears jewelry or watches under gloves?			□ Yes (-1)	□ No (+1)	Pts:		
Knows and demonstrates correct way to power-up BSC?			□ Yes (+1)	□ No (-1)	Pts:		
Resolves improper window slash height (or alarm)?			□ Yes (+1)	□ No (-1)	Pts:		
Identifies the	he BSC air-speed safe rar	□ Yes (+1)	□ No (-1)	Pts:			
Decontami	nates inside of BSC before	□ Yes (+2)	□ No (-2)	Pts:			
Establishes	s three essential work-area	□ Yes (+1)	□ No (-1)	Pts:			
Cleans hands with ethanol prior to introducing into BSC?			□ Yes (+1)	□ No (-1)	Pts:		
Cleans material with ethanol prior to introducing into BSC?			□ Yes (+1)	□ No (-1)	Pts:		
Introduces ethanol, NaOCl and waste bag inside BSC?			□ Yes (+1)	□ No (-1)	Pts:		
Introduces micropipettes, pipette tips and tubes inside BSC?			□ Yes (+1)	□ No (-1)	Pts:		
Places absorbent towel underneath NaOCl bottle?			□ Yes (+1)	□ No (-1)	Pts:		
Decontaminates inside of BSC before work?			□ Yes (+2)	□ No (-2)	Pts:		
Provides time for laminar flow to stabilize before work?			□ Yes (+1)	□ No (-1)	Pts:		





Adjusts chair for proper armpit height?	□ Yes (+1)	□ No (-1)	Pts:		
Movements inside BSC are slow and intentional?	□ Yes (+1)	□ No (-1)	Pts:		
Hands remain inside BSC at all times?	□ Yes (+1)	□ No (-1)	Pts:		
Front ventilations grills remain unobstructed at all times?	□ Yes (+1)	□ No (-1)	Pts:		
Decontaminates gloves before removing from BSC?	□ Yes (+1)	□ No (-1)	Pts:		
Removes outer gloves inside BSC after work in BSC?	□ Yes (+1)	□ No (-1)	Pts:		
Decontaminates material before removing from BSC?	□ Yes (+1)	□ No (-1)	Pts:		
Closes waste bag before removing from BSC?	□ Yes (+1)	□ No (-1)	Pts:		
Has to abandon BSC work due to forgetfulness?	□ Yes (+1)	□ No (-1)	Pts:		
Incurs in spills while working in BSC?	□ Yes (-2)	□ No (+2)	Pts:		
Removes NaOCl bottle after work in BSC?	□ Yes (+1)	□ No (-1)	Pts:		
Decontaminates inside of BSC after work?	□ Yes (+2)	□ No (-2)	Pts:		
Decontaminates inner gloves outside BSC after work in BSC?	□ Yes (+1)	□ No (-1)	Pts:		
Leaves UV light on after work in BSC?	□ Yes (+1)	□ No (-1)	Pts:		
Demonstrates correct way to shut-down BSC?	□ Yes (+1)	□ No (-1)	Pts:		
Describes Cell Biology lab decontamination procedure?	□ Yes (+1)	□ No (-1)	Pts:		
		Total Pt	s:/ 42		
Assessment result: □ Failed (≤29) □ Passed (≥30, 76	0%) □ Certified (≥36, 85%)				
ertified until:(dd/mmm/yyyy) NOTE: 1 year for those Passed, 3 years for those Certified					
Date of next assessment: (dd/mmm/yyyy)					

*Dr. Christian A. García-Sepúlveda*Principal Investigator
Viral & Human Genomics Lab
Facultad de Medicina UASLP

*Dr. Sandra E. Guerra-Palomares*Associate Investigator
Viral & Human Genomics Lab
Facultad de Medicina UASLP

