



## RT-PCR amplification of mammalian 18S rRNA

Created: Apr 03, 2019; Last modified: Jul 23, 2021 Version: 2.0

This protocol describes the procedure for 18S rRNA amplification for a wide array of molecular applications when working with eukaryotes. 18 Svedberg ribosomal RNA (18S rRNA) is the only rRNA molecule making up the small eukaryotic ribosomal subunit. 18S rRNA is a component of the small eukaryotic ribosomal subunit (40S) and a homologue of prokaryotes and mitochondrial 16S rRNA. The 18S rRNA genes are widely used in molecular analysis to reconstruct phylogeny but also made their way into molecular methods as internal control templates for both PCR and RT-PCR. This protocol has been specifically optimized for evaluating the PCR functionality of extracted RNA following other quality evaluations such as spectrophotometric and electrophoretic integrity assays. This protocol can be easily translated for work with genomic DNA for internal PCR controls and the like.

### Oligonucleotide primers

Name	Sequence*	bp	%GC	Tm	Hair	HmD	HtD	Amplicon	Ref
18S-F	5'-CgA-CgA-CCC-ATT-CgA-ACg-TCT-3'	21	57.1	59.3	-2.61	-10.65	-5.02	312 bp	1
18S-R	5'-gCT-ATT-ggA-gCT-ggA-ATT-ACC-g-3'	22	50	55.7	-1.35	-6.34			

\* Reverse oligonucleotide primer sequences given in this table are the reverse-complement of the sequence present in the alignments (see below) and as they should be ordered for synthesis. Hairpin (Hair), homodimer (Hm) and heterodimer (Ht)  $\Delta Q$ 's are given in kcal/mol.

### Oligonucleotide specificity

```

Rattus_norvegicus      CGACGACCCATTCGAACGTC-----GCCCTATCAACTTTCGATGGTAGTCGCCGTGCCTACCATGGT
Mus_musculus          -----
Homo_sapiens           -----
Drosophila_melanogaster -----AGAT T A T-----T-----AT-TAG-A-----
Aedes_aegypti         -----AGAT T A TA-----A-T-----ATAGAG-A-----
Aedes_albopictus      -----AGAT T A TA-----A-T-----ATAGAG-A-----
Culex_pipiens         -----AGAT T A TA-----A-T-----ATAGAG-A-----
Anopheles_pseudopunctipennis -----AGAT TCG A T-----A-T-----AT-GAG-A-----
Musa_coccinea         -----CAT A TT-----GATAG-G-----

Rattus_norvegicus      GACCACGGGTGACGGGGAATCAGGGTTCGATTCCGGAGAGGGAGCCTGAGAAACGGCTACCAC
Mus_musculus          -----
Homo_sapiens           -----
Drosophila_melanogaster TG-A-----A-----
Aedes_aegypti         TG-A-----A-----T-----
Aedes_albopictus      TG-A-----A-----T-----
Culex_pipiens         TG-A-----A-----T-----
Anopheles_pseudopunctipennis T--A-----A-----T-----
Musa_coccinea         -GTG-----A-----T-----

```





```

Rattus_norvegicus      ATCCAAGGAAGGCAGCAGGCGCGCAAATTACCCACTCCCGACCCGGGGAGGTAGTGACGAAAA
Mus_musculus          -----
Homo_sapiens          -----
Drosophila_melanogaster ---T-----T-----AG-T-----
Aedes_aegypti         -----T-----A---G-A-----G--
Aedes_albopictus      -----T-----A---G-A-----G--
Culex_pipiens         -----T-----A---G-A-----G--
Anopheles_pseudopunctipennis -----T-----A---G-A-----G--
Musa_coccinea         -----A--T---A-----A-T--

Rattus_norvegicus      ATAACAATACAGGACTCTTT••CGAGGCCCTGTAATTGGAATGAGTCCACTTTAAATCCTTTA
Mus_musculus          -----
Homo_sapiens          -----
Drosophila_melanogaster -----A-ATC-----A-----
Aedes_aegypti         -----T-A-----TAT--C-T-T-A-----TG-GCA-----C-
Aedes_albopictus      -----T-A-----TAT--C-T-T-A-----TG-GCA-----C-
Culex_pipiens         -----T-A-----TAT--T-TTT-A-----AC-G-GCA-----CG
Anopheles_pseudopunctipennis -----TGAA-----AAT--T-TTTCA-----TG-GCA-----
Musa_coccinea         -----C--G-----C--••GA-T-TG-----A--A-C-----C--

Rattus_norvegicus      ACGAGGATCCATTGGAGGGCAAGTCTGGTGCCAGCAGCCGCGGTAATTCAGCTCCAATAGC
Mus_musculus          -----
Homo_sapiens          -----
Drosophila_melanogaster --A--C-A-----
Aedes_aegypti         --A-----A-G-----C
Aedes_albopictus      --A-----A-G-----C
Culex_pipiens         GTA-----A-G-----C
Anopheles_pseudopunctipennis G-A-----A-G-----C
Musa_coccinea         -----
  
```

SURe v1.0 Sequence Unanimity Reformatting by UASLP's Viral & Human Genomics Laboratory (<http://midasmap.uaslp.mx/suretool/>) of Clustal Omega (1.2.4) multiple sequence alignment of GenBank retrieved accessions for *Mus musculus* (X00686.1), *Rattus norvegicus* (X01117.1), *Homo sapiens* (M10098.1), *Drosophila melanogaster* (M21017.1), *Aedes aegypti* (HE613439.1), *Aedes albopictus* (AB085210.1), *Culex pipiens* (HG001282.1), *Anopheles pseudopunctipennis* (U49735.1) and *Musa coccinea* (LC610755.1).

### Dual (separate RT and PCR) components and conditions



Distributed through a Creative Commons Attribution (BY) license granting the licensee the right to copy, distribute, display and make derivative works based on this document, including commercial use, as long as they credit the author as "Laboratorio de Genómica Viral y Humana, Facultad de Medicina UASLP".



Dual (separate RT and PCR) uses Vivantis M-MULV (Cat. No.: ME2305) RT Enzyme and Vivantis Taq DNA Polymerase (Cat. No.: PL1202), Vivantis Technologies Sdn Bhd, Malaysia

### First strand synthesis (RT)

		1x
dH <sub>2</sub> O	cf	2.6 µL
5 µM Forward oligo	1.125 µM	2.25 µL
5 µM Reverse oligo	1.125 µM	2.25 µL
10 mM dNTPs	250 µM	0.5 µL
RNA	-	10 µL
		vf: 17.6 µl

Total time: 6 min		
95 °C	2 min	1 cycle
4 °C	2 min	

Run RT-1 program in Axygen TC-1

RT Buffer	1x cf	2 µL
RT Enzyme 250 IU/µl	5 IU/µL	0.4 µL
		vf: 20 µl

Total time: 1:12 hrs		
38 °C	60 min	1 cycle
95 °C	5 min	
4 °C	5 min	

Run RT-2 program in Axygen TC-1

### Polymerase Chain Reaction (PCR)

		1x
dH <sub>2</sub> O	cf	7.58 µL
2x Buffer	1 x	1 µL
10 µM Forward oligo	150 nM	0.15 µL
50 mM MgCl <sub>2</sub>	1.5 mM	0.3 µL
10 mM dNTPs	200 µM	0.2 µL
5 µM oligo mix	200 nM	0.4 µL
5IU Taq DNA Pol	0.6 IU	0.12 µL
100 ng/µL cDNA	10 ng	0.4 µL
		vf: 10 µl


Total time: 2:10 hrs		
94 °C	2 min	35 cycles
94 °C	45 sec	
58 °C	30 sec	
72 °C	90 sec	
72 °C	10 min	

Run generic program in Applied Biosystems 7500



## Electrophoresis conditions

Load 12  $\mu\text{L}$  of 5  $\mu\text{L}$  of orange loading buffer plus the complete PCR reaction volume (10  $\mu\text{L}$ ), mix by pipetting, use 2.5% agarose gel in TAE buffer and run at 80 VDC for 75 minutes.



First lane corresponds to 100 bp molecular weight marker (MWM), NC corresponds to a  $\text{dH}_2\text{O}$  negative control, the following eight lanes correspond to duplicates of BALB/c intestinal tissue RNA extract. NOTE: Band present in NC corresponds to well 2 (BALB/c-01) overflow during gel loading.

## Notes

1. All RT-PCR reactions and mastermixes should be set up in the RT-PCR area ONLY!
2. Clean workbench with 0.1% NaOCl and 70% Ethanol before and after setting up the RT-PCR.
3. Preparation of all mastermixes should only be performed in the RT-PCR workbench.
4. Addition of sample RNA/DNA should only be performed in the RT-PCR workbench.
5. Prepare RT-PCR reactions on ice to prevent evaporation of small liquid volumes.
6. Vortex all reagents (except RNA, DNA and RT/Taq DNA polymerase) before dispensing into master mix.
7. Vortex master mix after adding all required reagents then distribute to each individual PCR tube.
8. Vortex all PCR reactions after adding reagents and DNA. Spin briefly and load into cycler.



## References

1. Zhou J. Microarrays for bacterial detection and microbial community analysis. *Curr Opin Microbiol.* 2003 Jun;6(3):288-94. Review. PMID: 12831906 [https://doi.org/10.1016/S1369-5274\(03\)00052-3](https://doi.org/10.1016/S1369-5274(03)00052-3).

## Revision history

- 1.0 Original document.
- 2.0 Expanded number of organisms included in alignments, changed title to indicate mammalian specificity of primers.

