



Middle East Respiratory Syndrome Coronavirus (MERS-CoV) detection & quantitation by qRT-PCR (SYBR green).

Created: May 23, 2017; Last modified: Nov 07, 2018; Version: 3.0

This protocol describes the real-time polymerase chain reaction with reverse transcription (qRT-PCR) approach developed for the detection and quantitation of Middle East Respiratory Syndrome Coronavirus (MERS-CoV) viral titres using a SYBR-Green modality (but easily TaqMan adaptable). CoV are a group IV ((+)ssRNA) species of the Nidovirales order, Coronaviridae family, Coronavirinae subfamily and Betacoronavirus genus having a twin species (SARS-CoV). MERS-CoV is a new member of Betacoronavirus, lineage C. MERS-CoV is distinct from SARS CoV and distinct from the common-cold CoV and known endemic human HCoV-OC43 and HCoV-HKU1. MERS-CoV is more closely related (more than 90% sequence identity) to bat coronaviruses HKU4 and HKU5 (lineage 2C) than to SARS-CoV (lineage 2B). MERS is considered to belong to the same species as bat CoV by the International Committee on Taxonomy of Viruses. The 26-32 kb CoV genome is non-segmented and contains 5'-methylated caps and 3'-polyadenylated tails encoding for replicase, structural proteins like spike glycoprotein (S), envelope protein (E), membrane protein (M) and nucleocapsid protein (N) which are essential for virus-cell-receptor binding and virion assembly, and immunomodulatory effects. The partially overlapping 5'-terminal ORF (ORF1a/b) encodes the large replicase polyprotein 1a (pp1a) and pp1ab which are cleaved by proteases to produce non-structural proteins, RNA-dependent RNA polymerase (RdRp) and helicase (Hel), involved in the transcription and replication of CoVs.

Oligonucleotide primers

Name	Size	bp	%GC	Tm	Hair	HomD	HetD	Size
MERS-F	5'-gCA-ACg-CgC-gAT-TCA-gTT-3'	18	55.6%	57°C	0.12	-17.11	-3.61	91
MERS-R	5'-gCC-TCT-ACA-Cgg-gAC-CCA-TA-3'	20	60%	58.5°C	-0.24	-6.14		

Note 1: Reverse oligonucleotide primer sequences given in this table are the reverse-complement of sequence present in alignments and as they should be ordered for synthesis; Hairpin (Hair), homodimer (homD) and heterodimer (HetD) ΔQ 's given in kcal/mol.

Oligonucleotide primer map in artificial gene² (underlined primer sites)

AAGCTCCGCGTCCTTTACCAAGGAAATAACATCACTGTAAGTGCCTATGCAAACGGCGACCATGCCGTCACAGTTAAGGACGCCA
AATTCATTGTGGGGCCAATGTCTTCAGCCTGGACACCTTTTCGACAACAAAATTGTGGTGTACAAAAGGTGACGTCTATAATATGGA
CTACCCGCCCTTTGGCGCAGGAAGACCAGGACAATTTGG | CAGACCACGCTACGGCGTGCTACTCTGCGGAGAGTGCAGTCTGCG
ATAGTGCCCGAGGAGACTGGGTTAACAAAGGCAAACCAACGC CCCACGCGGCCCTAG | AAGGACTAGAGGTTAGAGGAGACCAT
GTCCATGTCACCCACGGTCATCATCTGGCATGTCTTGGGTCTATCCTGCTGCTCTACAGCATCATTCAGGCACAGAACGCC
| AAGGACTAGAGGTTAGAGGAGACCCCCCAACACAAAACAGCATATTTGACGCTGGGAAAAGACCAGAGATCCTGCTGTCTCTGC
AACATCAATCCAGGCACAGAGCGCC | AAATACACATACCAAAAACAAAGTGGTGAAGGTTCTCAGACCAGCTGAAGGAGGGAAAAC
AGTCATGGACATCATCTCAAGACAAGACCAGAGAGGGAGTGGG | ATGTATGTGAGTGTGATGCCACGAAATGGTCAACAGGAGA
TAAATCGGCAAAGTTAAGAGATTCACACAGGCATTATATGATGGCTTGTGATGAGAA GTTAAAAATGTTGCGTTGTTGATGC |
ATGTAACACCTCTACAATG GATG CCGACAAGATTGTATTCAAAGTCAATAATCAGGTGGTCTCTTTGAAGCCTGAGATTATCGTG
GATCAAATGAGTACAAGTACCCCTGC | GCATAGCTGTAGGAAGGACTCTCAAAGGCATTTTTCTCTGCTGCTTACTGACCCCTTT





AGGGAACGAAGCCCCTGGGGGATACTGTCTTGAAAAGTGGATGCTT | GCAACGCGCGATTTCAGTTCTCTTCACATAATCGCCCC
GAGCTCGCTTATCGTTTAAAGCAGCTCTGCGCTACTATGGGTCCCCTGTAGAGGC | TGATGATGCCGTCGTGTGCTACAA CAGTAA
CTATGCGGCTCAAGGTTTAGTAGCTAGCATTAAAGAACTTAAAGGCAGTTCTGTATTATCAAGATAATGTGTTTCATGTCTGAGGCA
AAATGTTGGACTGAGACTGACCTTACTAAA GGACCTCACGAATTTTGCTCACA

Note 2: Primer sites for different pathogen sequences shown in red, TaqMan probe sites shown in blue, MERS-CoV specific primer binding sites shown in grey highlight. MERS-CoV artificial gene sequence derived from upE ORF 5 gene sequence as suggested by Corman 2012.

Two-step (separate RT and PCR) components and conditions

Dual (separate RT and PCR) uses M-MULV RT Enzyme (Vivantis Cat. No.: ME2305) and iTaq qPCR Master Mix (Biorad Cat. No.: 1725150).

First strand synthesis (RT)

		1x
dH ₂ O	cf	2.6 µL
10 µM Forward oligo	1.125 µM	2.25 µL
10 µM Reverse oligo	1.125 µM	2.25 µL
10 mM dNTPs 10 mM	250 µM	0.5 µL
RNA	-	10 µL
		vf: 17.6 µl

Total time: 6 min		
95 °C	2 min	1 cycle
4 °C	2 min	

Run RT-1 program in Axygen TC-1

RT Buffer	1x cf	2 µL
RT Enzyme 250 IU/µl	5 IU/µL	0.4 µL
		vf: 20 µl

Total time: 1:12 hrs		
38 °C	60 min	1 cycle
95 °C	5 min	
4 °C	5 min	

Run RT-2 program in Axygen TC-1

Polymerase Chain Reaction (PCR)

		1x
dH ₂ O	cf	3.7 µL
Master mix	1 x	5 µL
10 µM Forward oligo	150 nM	0.15 µL
10 µM Reverse oligo	150 nM	0.15 µL
Template	10 ng	1 µL
		vf: 10 µl

Total time: 2:10 hrs		
95 °C	3 min	40 cycles
94 °C	15 sec	
60 °C ³	60 sec	
95 °C	15 sec	
60 °C	20 sec	
Ramp ³	60 min	
95 °C	15 sec	

Run generic program in Applied Biosystems 7500

Note 3: Data acquisition.



One-step qRT-PCR components and conditions

The One-Step qRT-PCR procedure uses VeriQuest SYBR Green One-Step qRT-PCR Mix (USB Affymetrix Cat. No.: 75700).

qRT-PCR		1x
dH ₂ O	cf	3.1 μL
qRT-PCR Master mix	1 x	5 μL
RT Enzyme mix	1 x	0.1 μL
10 μM Forward oligo	400 nM	0.4 μL
10 μM Reverse oligo	400 nM	0.4 μL
Template	< 1 μg	1 μL
		vf: 10 μl

↓

Run Affymetrix program in Applied 7500		
--	--	--

→

Total time: 2:00 hrs		
50 °C	10 min	40 cycles
95 °C	10 min	
95 °C	15 sec	
58 °C ³	30 sec	
95 °C	15 sec	
60 °C	20 sec	
Ramp ³	60 min	
95 °C	15 sec	

Note 3: Data acquisition.



Titration curve preparation for quantitative analysis of viral titres

Using Artificial Gene (AG) sequence stock having all viral references (CHIKV, WNV, DENV1-3, DENV2-4, ZIKV, SNV, RABV, TCRV, MERS and SARS) prepare the following log titres (stock currently at 7.51×10^7 cp/μL, working dilution of AG prepared at 7.51×10^6 cp/μL).

	AG stock	dH ₂ O vol	Ct mean	
1 st	7.5×10^7 cp/μL	10 μL	90 μL	7.70 (±0.74)
2 nd	7.5×10^6 cp/μL	10 μL	90 μL	12.63 (±0.34)
3 rd	7.5×10^5 cp/μL	10 μL	90 μL	16.48 (±0.17)
4 th	7.5×10^4 cp/μL	10 μL	90 μL	19.52 (±0.24)
5 th	7.5×10^3 cp/μL	10 μL	90 μL	23.82 (±0.72)
6 th	7.5×10^2 cp/μL	10 μL	90 μL	26.64 (±0.08)

Add 90 μL to each of the 6 PCR 0.2 mL tubes. Take 10 μL of initial working stock (at 7.51×10^7 cp/μL) and dispense into 1st PCR tube, wash tip 30 times, cap, vortex for 10 seconds and spin down for 10 seconds. Retrieve 10 μL from volumetric centre of PCR tube and dispense into 2nd PCR tube repeating exactly the same procedure for further dilutions.

Performance summary



Standard curve: $m = -3.75$, $Y = 37.955$, $R^2 = 0.99$
 Target amplicon T_m : $80.18^\circ\text{C} \pm 0.17 \text{ SD}$ (80.01 to 80.35°C)
 Limit of detection (LODet): $7.51 \times 10^0 \text{ cp}/\mu\text{L}$
 Limit of discrimination (LODis): $7.51 \times 10^2 \text{ cp}/\mu\text{L}$

Interpretation

For samples to be considered as positive, their C_t should not be prior to that of the last standard curve's C_t (i.e., $< \text{cycle } 27$, corresponding to the 7.5×10^2 LODet titre. In addition, the melting curve of the corresponding sample should be higher than surrounding unspecific melting peaks and preferably clean (without accompanying shoulders or unspecific melting peaks).

References

1. Corman VM, Eckerle I, Bleicker T, Zaki A, Landt O, Eschbach-Bludau M, van Boheemen S, Gopal R, Ballhause M, Bestebroer TM, Muth D, Müller MA, Drexler JF, Zambon M, Osterhaus AD, Fouchier RM, Drosten C. Detection of a novel human coronavirus by real-time reverse-transcription polymerase chain reaction. *Euro Surveill.* 2012;17(39).
2. Zumla A, Chan JF, Azhar EI, Hui DS, Yuen KY. Coronaviruses - drug discovery and therapeutic options. *Nat Rev Drug Discov.* 2016 May;15(5):327-47. doi: 10.1038/nrd.2015.37. Epub 2016 Feb 12. Review. PubMed PMID: 26868298.



Revision history

- 1.0 Original document.
- 2.0 Optimized protocol.
- 3.0 Changes to document format only.

